




EFCB7 rabbit pAb

Catalog No	YP-Ab-08253
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	EFCAB7 KIAA1799
Immunogen	Synthesized peptide derived from human EFCB7 AA range: 386-436
Specificity	This antibody detects endogenous levels of EFCB7 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EFCB7;EFCAB7 KIAA1799
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell projection, cilium membrane ; Peripheral membrane protein ; Cytoplasmic side . The EvC complex localizes at the base of cilia in the EvC zone of primary cilia in a EFCAB7-dependent manner. .
Tissue Specificity	
Function	similarity:Contains 3 EF-hand domains.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from K562 cells,
primary antibody was diluted at 1:1000, 4° over night