




CI163 rabbit pAb

Catalog No	YP-Ab-08257
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	C9orf163
Immunogen	Synthesized peptide derived from human CI163 AA range: 144-194
Specificity	This antibody detects endogenous levels of CI163 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CI163;C9orf163
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	miscellaneous:This protein is produced by a bicistronic transcript which also produces the C9orf163 protein from an non-overlapping reading frame., PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from 293T cells,
primary antibody was diluted at 1:1000, 4° over night