




ZN277 rabbit pAb

Catalog No	YP-Ab-08262
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	ZNF277 NRIF4 ZNF277P
Immunogen	Synthesized peptide derived from human ZN277 AA range: 201-251
Specificity	This antibody detects endogenous levels of ZN277 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NRIF 4;NRIF4;Nuclear receptor interacting factor 4;Nuclear receptor-interacting factor 4;Zinc finger protein (C2H2 type) 277;Zinc finger protein 277;Zinc finger protein 277 pseudogene;ZN277_HUMAN;ZNF 277;ZNF277;ZNF277P;ZN277; ZNF277 NRIF4 ZNF277P
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	
Function	function:May be involved in transcriptional regulation.,similarity:Belongs to the ZNF277 family.,similarity:Contains 2 C2H2-type zinc fingers.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from A549 cells,
primary antibody was diluted at 1:1000, 4° over night