




AKND1 rabbit pAb

Catalog No	YP-Ab-08266
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	AKNAD1 C1orf62
Immunogen	Synthesized peptide derived from human AKND1 AA range: 373-423
Specificity	This antibody detects endogenous levels of AKND1 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AKND1;AKNAD1 C1orf62
Observed Band	
Calculated Molecular Weight	
Cell Pathway	cytoplasm,
Tissue Specificity	
Function	
Background	This gene encodes a protein which contains a domain found in an AT-hook-containing transcription factor. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2012],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HuvEc cells,
primary antibody was diluted at 1:1000, 4° over night