




GPHRB rabbit pAb

Catalog No	YP-Ab-08271
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	GPR89B GPHRB HSPC201
Immunogen	Synthesized peptide derived from human GPHRB AA range: 270-320
Specificity	This antibody detects endogenous levels of GPHRB at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GPHRB;GPR89B GPHRB HSPC201
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Golgi apparatus membrane ; Multi-pass membrane protein .
Tissue Specificity	Ubiquitous.
Function	caution:Could be the product of a pseudogene.,function:May be involved in the NF-kappa-B and MAPK signaling pathways.,miscellaneous:Does not seem to be able to bind GTP.,sequence caution:Translated as Cys.,similarity:Belongs to the GPR89 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from SH-SY5Y cells,
primary antibody was diluted at 1:1000, 4°over night