




ARMC1 rabbit pAb

Catalog No	YP-Ab-08272
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	ARMC1 ARCP
Immunogen	Synthesized peptide derived from human ARMC1 AA range: 203-253
Specificity	This antibody detects endogenous levels of ARMC1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SVH;PNAS112;PNAS-112;PSEC0198;ARMC10;ARMC1;ARMC1 ARCP
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Mitochondrion . Mitochondrion outer membrane . Associates with the outer mitochondrion membrane, most likely through its C-terminus (PubMed:31644573). Not integrated into the mitochondrial outer membrane (PubMed:31644573). .
Tissue Specificity	
Function	similarity:Contains 1 ARM repeat.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night