




DYRL1 rabbit pAb

Catalog No	YP-Ab-08286
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	DHFRL1 DHFRP4
Immunogen	Synthesized peptide derived from human DYRL1 AA range: 31-81
Specificity	This antibody detects endogenous levels of DYRL1 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DYRL1;DHFRL1 DHFRP4
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Mitochondrion . Mitochondrion matrix . Mitochondrion inner membrane .
Tissue Specificity	Expressed in numerous cell lines.
Function	similarity:Belongs to the dihydrofolate reductase family.,similarity:Contains 1 DHFR (dihydrofolate reductase) domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from PC-12 cells,
primary antibody was diluted at 1:1000, 4° over night