




SCND3 rabbit pAb

Catalog No	YP-Ab-08295
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	SCAND3 KIAA1925 ZNF305P2 ZNF452
Immunogen	Synthesized peptide derived from human SCND3 AA range: 1111-1161
Specificity	This antibody detects endogenous levels of SCND3 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SCND3;SCAND3 KIAA1925 ZNF305P2 ZNF452
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	Weakly expressed in the lung (at protein level).
Function	similarity:Contains 1 integrase catalytic domain.,similarity:Contains 1 SCAN box domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from SW480 cells,
primary antibody was diluted at 1:1000, 4° over night