




GINM1 rabbit pAb

Catalog No	YP-Ab-08297
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	GINM1 C6orf72 UNQ710/PRO1361
Immunogen	Synthesized peptide derived from human GINM1 AA range: 241-291
Specificity	This antibody detects endogenous levels of GINM1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GINM1;GINM1 C6orf72 UNQ710/PRO1361
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Single-pass type I membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from SH-SY5Y cells,
primary antibody was diluted at 1:1000, 4° over night