



CAN15 rabbit pAb


Catalog No	YP-Ab-08304
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	SOLH CAPN15
Immunogen	Synthesized peptide derived from human CAN15 AA range: 153-203
Specificity	This antibody detects endogenous levels of CAN15 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	281411K16Rik;MGC1325;Zinc finger and SCAN domain containing protein 15; Zinc finger and SCAN domain-containing protein 15;Zinc finger protein 397;Zinc finger protein 47;ZN397_HUMAN;ZNF 397;ZNF 47;ZNF397;ZNF47;ZSCAN 15; ZSCAN15;CAN15;SOLH CAPN15
Observed Band	
Calculated Molecular Weight	
Cell Pathway	cytoplasm,
Tissue Specificity	Widely expressed with higher expression in brain.
Function	similarity:Belongs to the peptidase C2 family.,similarity:Contains 1 calpain catalytic domain.,similarity:Contains 5 RanBP2-type zinc fingers.,tissue specificity:Widely expressed with higher expression in brain.,
Background	This gene encodes a protein containing zinc-finger-like repeats and a calpain-like protease domain. The encoded protein may function as a transcription factor, RNA-binding protein, or in protein-protein interactions during visual system development. [provided by RefSeq, Jul 2008],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products ImagesWestern blot analysis of lysates from HpeG2 cells,
primary antibody was diluted at 1:1000, 4° over night