




BIRC8 rabbit pAb

Catalog No	YP-Ab-08307
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	BIRC8 ILP2
Immunogen	Synthesized peptide derived from human BIRC8 AA range: 122-172
Specificity	This antibody detects endogenous levels of BIRC8 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Baculoviral IAP repeat containing 8;Baculoviral IAP repeat containing protein 8; Baculoviral IAP repeat-containing protein 8;BIRC 8;BIRC8;BIRC8 protein; BIRC8_HUMAN;hILP 2;hILP2;IAP like protein 2;IAP-like protein 2;ILP 2;ILP-2; Inhibitor of apoptosis like protein 2;Inhibitor of apoptosis-like protein 2;Testis specific inhibitor of apoptosis;Testis-specific inhibitor of apoptosis;BIRC8 ILP2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm .
Tissue Specificity	Testis specific in normal tissues.
Function	function:Protects against apoptosis mediated by BAX.,similarity:Belongs to the IAP family.,similarity:Contains 1 BIR repeat.,similarity:Contains 1 RING-type zinc finger.,subunit:Binds to caspase-9.,tissue specificity:Testis specific in normal tissues.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from CACO2 cells, primary antibody was diluted at 1:1000, 4° over night