




PTH rabbit pAb

Catalog No	YP-Ab-08309
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	PTRH1 C9orf115
Immunogen	Synthesized peptide derived from human PTH AA range: 53-103
Specificity	This antibody detects endogenous levels of PTH at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	parathyroid hormone;PTH;Parathormone;Parathyrin;PTRH1 C9orf115
Observed Band	
Calculated Molecular Weight	
Cell Pathway	mitochondrion,
Tissue Specificity	
Function	catalytic activity:N-substituted aminoacyl-tRNA + H(2)O = N-substituted amino acid + tRNA.,similarity:Belongs to the PTH family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from Hela cells,
primary antibody was diluted at 1:1000, 4° over night