




ANM10 rabbit pAb

Catalog No	YP-Ab-08315
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	PRMT10
Immunogen	Synthesized peptide derived from human ANM10 AA range: 366-416
Specificity	This antibody detects endogenous levels of ANM10 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ANM1_HUMAN;Prmt1;protein arginine methyltransferase 1 (putative);Putative protein arginine N-methyltransferase 1;ANM10;PRMT10
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from A431 cells,
primary antibody was diluted at 1:1000, 4° over night