




CJ126 rabbit pAb

Catalog No	YP-Ab-08316
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	C10orf126
Immunogen	Synthesized peptide derived from human CJ126 AA range: 78-128
Specificity	This antibody detects endogenous levels of CJ126 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CJ126;C10orf126
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	sequence caution:Wrong choice of frame.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night