




T132D rabbit pAb

Catalog No	YP-Ab-08325
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	TMEM132D HBE120 KIAA1944 MOLT
Immunogen	Synthesized peptide derived from human T132D AA range: 509-559
Specificity	This antibody detects endogenous levels of T132D at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	T132D;TMEM132D HBE120 KIAA1944 MOLT
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Single-pass type I membrane protein .
Tissue Specificity	Expressed in mature oligodendrocytes. Detected in the brain, lung, pancreas and testis.
Function	function:May serve as a cell-surface marker for oligodendrocyte differentiation., similarity:Belongs to the TMEM132 family.,tissue specificity:Expressed in mature oligodendrocytes. Detected in the brain, lung, pancreas and testis.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from AD293 cells,
primary antibody was diluted at 1:1000, 4°over night