



T132E rabbit pAb


Catalog No	YP-Ab-08331
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	TMEM132E
Immunogen	Synthesized peptide derived from human T132E AA range: 486-536
Specificity	This antibody detects endogenous levels of T132E at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	T132E;TMEM132E
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Single-pass type I membrane protein .
Tissue Specificity	
Function	disease:TMEM132E gene is identified as one gene located in a region involved in a heterozygous deletion of approximately 4.7 Mb; this deletion, involving the NF1 gene and contiguous genes lying in its flanking regions, is observed in one patient with NF1 and 17q11.2 microdeletion syndrome [MIM:162200]. The NF1 microdeletion syndrome is often characterized by a more severe phenotype than that observed in the majority of NF1 patients. Indeed, patients with NF1 microdeletion often show variable facial dysmorphism, mental retardation, developmental delay, and an excessive number of neurofibromas., similarity:Belongs to the TMEM132 family.,
Background	

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from HEK293 cells, primary antibody was diluted at 1:1000, 4° over night