



CD11A rabbit pAb


Catalog No	YP-Ab-08340
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	CDK11A CDC2L2 CDC2L3 PITSLREB
Immunogen	Synthesized peptide derived from human CD11A AA range: 53-103
Specificity	This antibody detects endogenous levels of CD11A at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CD11 antigen-like family member A;Integrin alpha L;Integrin alpha-L;ITGAL; Leukocyte adhesion glycoprotein LFA-1 alpha chain;CD11A;CDK11A CDC2L2 CDC2L3 PITSLREB
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm. Nucleus.
Tissue Specificity	Expressed ubiquitously. Some evidence of isoform-specific tissue distribution.
Function	caution:The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data.,
Background	This gene encodes a member of the serine/threonine protein kinase family. Members of this kinase family are known to be essential for eukaryotic cell cycle control. Due to a segmental duplication, this gene shares very high sequence identity with a neighboring gene. These two genes are frequently deleted or altered in neuroblastoma. The protein kinase encoded by this gene can be cleaved by caspases and may play a role in cell apoptosis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2015],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products ImagesWestern blot analysis of lysates from HeLa cells,
primary antibody was diluted at 1:1000, 4°over night