




LRRT2 rabbit pAb

Catalog No	YP-Ab-08343
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	LRRTM2 KIAA0416 LRRN2
Immunogen	Synthesized peptide derived from human LRRT2 AA range: 372-422
Specificity	This antibody detects endogenous levels of LRRT2 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LRRT2; LRRTM2 KIAA0416 LRRN2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane; Single-pass type I membrane protein. Cell junction, synapse, postsynaptic cell membrane ; Single-pass type I membrane protein . Localized to excitatory synapses. .
Tissue Specificity	Expressed in neuronal tissues.
Function	function: May play a role in the development and maintenance of the vertebrate nervous system., similarity: Belongs to the LRRTM family., similarity: Contains 10 LRR (leucine-rich) repeats., tissue specificity: Expressed in neuronal tissues.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from 3T3 cells, primary antibody was diluted at 1:1000, 4° over night