




GET4 rabbit pAb

Catalog No	YP-Ab-08351
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	GET4 C7orf20 CEE TRC35 CGI-20
Immunogen	Synthesized peptide derived from human GET4 AA range: 190-240
Specificity	This antibody detects endogenous levels of GET4 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C7orf20;CEE;CGI 20;CGI-20;Conserved edge-expressed protein;GET4;GET4 C7orf20 CEE TRC35 CGI-20
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, cytosol .
Tissue Specificity	
Function	similarity:Belongs to the UPF0363 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from A431 cells,
primary antibody was diluted at 1:1000, 4° over night