




TM63B rabbit pAb

Catalog No	YP-Ab-08352
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	TMEM63B C6orf110
Immunogen	Synthesized peptide derived from human TM63B AA range: 627-677
Specificity	This antibody detects endogenous levels of TM63B at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TM63B;TMEM63B C6orf110
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	similarity:Belongs to the SPO75/TMEM63 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from 293T cells,
primary antibody was diluted at 1:1000, 4° over night