



USBP1 rabbit pAb


Catalog No	YP-Ab-08358
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	USHBP1 AIEBP MCC2
Immunogen	Synthesized peptide derived from human USBP1 AA range: 185-235
Specificity	This antibody detects endogenous levels of USBP1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AIE 75 binding protein; AIE 75 binding protein protein; AIE-75-binding protein; AIEBP; MCC 2; MCC-2; MCC2; Mutated in colon cancer protein 2; USBP1_HUMAN; USH1C binding protein 1; USH1C-binding protein 1; Ushbp1; Usher syndrome 1C binding protein 1; Usher syndrome type 1C protein binding protein 1; Usher syndrome type-1C protein-binding protein 1; USBP1; USHBP1 AIEBP MCC2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Highest level of expression in heart, and moderate to low expression in skeletal muscle, kidney, liver, small intestine, placenta and lung.
Function	similarity: Belongs to the MCC family., subunit: Interacts via its C-terminus with the first PDZ domain of USH1C., tissue specificity: Highest level of expression in heart, and moderate to low expression in skeletal muscle, kidney, liver, small intestine, placenta and lung.,
Background	

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from MDA-MB cells, primary antibody was diluted at 1:1000, 4° over night