




YN3516

CT173 rabbit pAb

Catalog No	YP-Ab-08359
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	C20orf173
Immunogen	Synthesized peptide derived from human CT173 AA range: 31-81
Specificity	This antibody detects endogenous levels of CT173 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CT173;C20orf173
Observed Band	
Calculated Molecular Weight	
Cell Pathway	integral component of membrane,
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from MCF-7 cells, primary antibody was diluted at 1:1000, 4° over night