




AKIP rabbit pAb

Catalog No	YP-Ab-08360
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	AURKAIP1 AIP AKIP
Immunogen	Synthesized peptide derived from human AKIP AA range: 10-60
Specificity	This antibody detects endogenous levels of AKIP at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	28S ribosomal protein S38;mitochondrial;AIP;AKIP;AURKA interacting protein; AURKA-interacting protein;AURKAIP1 AIP AKIP
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Mitochondrion . Nucleus .
Tissue Specificity	Ubiquitously expressed and especially highly expressed in heart, skeletal muscle and testis.
Function	function:May act as a negative regulator of Aurora-A kinase, by down-regulation through proteasome-dependent degradation.,subunit:Interacts with Aurora-A., tissue specificity:Ubiquitously expressed and especially highly expressed in heart, skeletal muscle and testis.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from THP-1 cells, primary antibody was diluted at 1:1000, 4° over night