




F162A rabbit pAb

Catalog No	YP-Ab-08364
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	FAM162A C3orf28 E2IG5 DC16 FWP001
Immunogen	Synthesized peptide derived from human F162A AA range: 78-128
Specificity	This antibody detects endogenous levels of F162A at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	G/G1 switch regulatory protein 24; GS24; GOS24; Growth factor-inducible nuclear protein NUP475; NUP475; Protein TIS11A; RNF162A; TIS 11; TIS11; TIS11A; Tristetraprolin; Tristetraproline; TTP; TTP_HUMAN; Zfp-36; ZFP36; Zinc finger protein 36; Zinc finger protein 36 homolog; C3H type; homolog (mouse); Zinc finger protein; 36 homolog; F162A; FAM162A C3orf28 E2IG5 DC16 FWP001
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Mitochondrion membrane ; Single-pass membrane protein .
Tissue Specificity	
Function	induction: By 17-beta-estradiol., similarity: Belongs to the UPF0389 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from HpeG2 cells, primary antibody was diluted at 1:1000, 4° over night