



DHR12 rabbit pAb

Catalog No	YP-Ab-08370
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	DHRS12
Immunogen	Synthesized peptide derived from human DHR12 AA range: 83-133
Specificity	This antibody detects endogenous levels of DHR12 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Alc alpha;Alc-alpha;Alcadein alpha 1;Alcadein alpha;Alcadein-alpha;alcalpha1;alcalpha2;Alzheimer related cadherin like protein;Alzheimer-related cadherin-like protein;C-terminal fragment 1-alpha;Cadherin related family member 12;Calsyntenin 1;Calsyntenin1;CDHR12;CLSTN 1;Clstn1;CS1;CSTN1;CSTN1_HUMAN;CTF1-alpha;FLJ32258;KIAA911;Non classical cadherin XB31alpha;Non classical cadherin XB31alpha1;Non-classical cadherin XB31alpha;PIK3CD;SAlc-alpha;XB31alpha;DHR12;DHRS12
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	function:Putative oxidoreductase.,similarity:Belongs to the short-chain dehydrogenases/reductases (SDR) family.,
Background	This gene encodes a member of the short-chain dehydrogenases/reductases (SDR) family, which has over 46,000 members. Members in this family are enzymes that metabolize many different compounds, such as steroid hormones, prostaglandins, retinoids, lipids and xenobiotics. Alternative splicing results in

multiple transcript variants and protein isoforms. [provided by RefSeq, Jul 2012],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from VEC cells,
primary antibody was diluted at 1:1000, 4° over night