




CRERF rabbit pAb

Catalog No	YP-Ab-08371
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	CREBRF C5orf41
Immunogen	Synthesized peptide derived from human CRERF AA range: 74-124
Specificity	This antibody detects endogenous levels of CRERF at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CRERF;CREBRF C5orf41
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus . Colocalizes with CREB3 in nuclear foci.
Tissue Specificity	
Function	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Belongs to the UPF0474 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MCF-7 cells,
primary antibody was diluted at 1:1000, 4° over night