




SNX7 rabbit pAb

Catalog No	YP-Ab-08376
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	SNX7
Immunogen	Synthesized peptide derived from human SNX7 AA range: 253-303
Specificity	This antibody detects endogenous levels of SNX7 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SNX7;sorting nexin 7
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Early endosome membrane ; Peripheral membrane protein ; Cytoplasmic side .
Tissue Specificity	
Function	function:May be involved in several stages of intracellular trafficking., similarity:Belongs to the sorting nexin family.,similarity:Contains 1 PX (phox homology) domain.,
Background	This gene encodes a member of the sorting nexin family. Members of this family contain a phox (PX) domain, which is a phosphoinositide binding domain, and are involved in intracellular trafficking. This protein does not contain a coiled coil region like some family members, and its exact function is unknown. Alternative splicing results in multiple transcript variants. A related pseudogene has been identified on chromosome 11. [provided by RefSeq, Jun 2010],
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from AD293 cells, primary antibody was diluted at 1:1000, 4° over night