




YN3536

CC080 rabbit pAb

Catalog No	YP-Ab-08379
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	C3orf80
Immunogen	Synthesized peptide derived from human CC080 AA range: 101-151
Specificity	This antibody detects endogenous levels of CC080 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CC080;C3orf80
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Single-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from HuvEc cells, primary antibody was diluted at 1:1000, 4° over night