



YN3539

LINES rabbit pAb

Catalog No	YP-Ab-08382
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	LINS LINS1 WINS1
Immunogen	Synthesized peptide derived from human LINES AA range: 105-155
Specificity	This antibody detects endogenous levels of LINES at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Acylcholine acylhydrolase;BCHE;Butyrylcholine esterase;CHE1;CHE2;CHLE_HUMAN;Choline esterase II;Cholinesterase (serum) 2;Cholinesterase 1;Cholinesterase;E1;Pseudocholinesterase;Pseudocholinesterase E1;LINES;LINS LINS1 WINS1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Expressed in adult testis, prostate, prostate, spleen, thymus, skeletal muscle, fetal kidney and brain.
Function	similarity:Belongs to the protein lines family.,tissue specificity:Expressed in adult testis, prostate, prostate, spleen, thymus, skeletal muscle, fetal kidney and brain.,



Background

The Drosophila segment polarity gene *lin* encodes a protein, *lines*, which plays important roles in development of the epidermis and hindgut. This gene encodes a protein containing a *lines*-like domain. This gene is located on chromosome 15 and clustered with the gene encoding ankyrin repeat and SOCS box-containing protein 7. [provided by RefSeq, Sep 2010],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from 293T cells, primary antibody was diluted at 1:1000, 4° over night