




CEP89 rabbit pAb

Catalog No	YP-Ab-08383
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	CEP89 CCDC123
Immunogen	Synthesized peptide derived from human CEP89 AA range: 238-288
Specificity	This antibody detects endogenous levels of CEP89 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CCDC123;CEP89;Centrosomal protein 123;Centrosomal protein of 89 kDa; Cep123;CEP89 CCDC123
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, cytosol. Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Cytoplasm, cytoskeleton, spindle pole. Cytoplasm, cytoskeleton, microtubule organizing center, centrosome, centriole. Mitochondrion intermembrane space. Localizes to the distal appendage region of the centriole, which anchors the mother centriole to the plasma membrane.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from HeLa cells, primary antibody was diluted at 1:1000, 4° over night