




YN3547

TM141 rabbit pAb

Catalog No	YP-Ab-08390
Isotype	IgG
Reactivity	Human; Mouse
Applications	IHC;IF
Gene Name	TMEM141
Immunogen	Synthesized peptide derived from human TM141 AA range: 5-55
Specificity	This antibody detects endogenous levels of TM141 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:50-200. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AHUS 6;AHUS6;BDCA 3;BDCA3;CD 141;CD141;CD141 antigen;Fetomodulin;Thbd;THPH12;THRM;Thrombomodulin;TM;TRBM_HUMAN;TM141;TMEM141
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	similarity:Belongs to the TMEM141 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Immunohistochemical analysis of paraffin-embedded human colon cancer. 1, Tris-EDTA, pH9.0 was used for antigen retrieval. 2 Antibody was diluted at 1:200 (4° overnight). 3, Secondary antibody was diluted at 1:200 (room temperature, 45min).