




TLX2 rabbit pAb

Catalog No	YP-Ab-08393
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	TLX2 HOX11L1 NCX
Immunogen	Synthesized peptide derived from human TLX2 AA range: 17-67
Specificity	This antibody detects endogenous levels of TLX2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Homeobox protein Hox-11L1;HOX11L1;NCX;T-cell leukemia homeobox protein 2;TLX 2;TLX2;TLX2 HOX11L1 NCX
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	
Function	similarity:Contains 1 homeobox DNA-binding domain.,
Background	This gene is a member of an orphan homeobox-containing transcription factor family. Studies of the mouse ortholog have shown that the encoded protein is crucial for the development of the enteric nervous system; in humans, loss-of-function may play a role in tumorigenesis of gastrointestinal stromal tumors. [provided by RefSeq, May 2010],
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from K562 cells,
primary antibody was diluted at 1:1000, 4°over night