




NAB1 rabbit pAb

Catalog No	YP-Ab-08394
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	NAB1
Immunogen	Synthesized peptide derived from human NAB1 AA range: 43-93
Specificity	This antibody detects endogenous levels of NAB1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EGR 1 binding protein 1; NAB1; NGFI A binding protein 1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	Isoform Short is found in myeloid leukemia cell line KG-1.
Function	domain: The NAB conserved domain 1 (NCD1) interacts with EGR1 inhibitory domain and mediates multimerization., domain: The NAB conserved domain 2 (NCD2) is necessary for transcriptional repression., function: Acts as a transcriptional repressor for zinc finger transcription factors EGR1 and EGR2., similarity: Belongs to the NAB family., subunit: Homomultimers may associate with EGR1 bound to DNA., tissue specificity: The short isoform is found in a myeloid Leukemia cell line.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from CACO2 cells, primary antibody was diluted at 1:1000, 4° over night