




GOGA7 rabbit pAb

Catalog No	YP-Ab-08395
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	GOLGA7 GCP16 HDCKB03P HSPC041
Immunogen	Synthesized peptide derived from human GOGA7 AA range: 64-114
Specificity	This antibody detects endogenous levels of GOGA7 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GOGA7; GOLGA7 GCP16 HDCKB03P HSPC041
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Golgi apparatus membrane ; Lipid-anchor .
Tissue Specificity	Expressed in all tissues except colon and thymus.
Function	function: May be involved in protein transport from Golgi to cell surface. The ZDHHC9-GOLGA7 complex is a palmitoyltransferase specific for HRAS and NRAS., PTM: Palmitoylated on Cys-69 and Cys-72; which is required for Golgi localization and interaction with GOLGA3., similarity: Belongs to the ERF4 family., subunit: Interacts with GOLGA3. Interacts with ZDHHC9., tissue specificity: Expressed in all tissues except colon and thymus.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from HpeG2 cells, primary antibody was diluted at 1:1000, 4° over night