




NFXL1 rabbit pAb

Catalog No	YP-Ab-08399
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	NFXL1 OZFP
Immunogen	Synthesized peptide derived from human NFXL1 AA range: 103-153
Specificity	This antibody detects endogenous levels of NFXL1 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	hOZFP;NF-X1-type zinc finger protein NFXL1;NFXL1;NFXL1_HUMAN;Nuclear transcription factor;Ovarian zinc finger protein;Up-regulated in colon cancer 5; URCC5;X-box binding-like 1;NFXL1 OZFP
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Single-pass membrane protein .
Tissue Specificity	
Function	similarity:Belongs to the NFX1 family.,similarity:Contains 1 RING-type zinc finger.,similarity:Contains 10 NF-X1-type zinc fingers.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from U2OS cells,
primary antibody was diluted at 1:1000, 4° over night