




MET22 rabbit pAb

Catalog No	YP-Ab-08403
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	METTL22 C16orf68 LP8272
Immunogen	Synthesized peptide derived from human MET22 AA range: 60-110
Specificity	This antibody detects endogenous levels of MET22 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MET22;METTL22 C16orf68 LP8272
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	
Function	
Background	This gene encodes a member of the non-histone lysine methyltransferases. It interacts with its substrate, Kin17, which is involved in DNA repair and replication and mRNA processing. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2016],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HCT116 cells,
primary antibody was diluted at 1:1000, 4°over night