




UAP1L rabbit pAb

Catalog No	YP-Ab-08404
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	UAP1L1
Immunogen	Synthesized peptide derived from human UAP1L AA range: 287-337
Specificity	This antibody detects endogenous levels of UAP1L at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	UAP1L;UAP1L1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	cytosol,
Tissue Specificity	
Function	similarity:Belongs to the UDGP type 1 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from AD293 cells,
primary antibody was diluted at 1:1000, 4° over night