



KTAP2 rabbit pAb


Catalog No	YP-Ab-08409
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	KRTCAP2 KCP2
Immunogen	Synthesized peptide derived from human KTAP2 AA range: 39-89
Specificity	This antibody detects endogenous levels of KTAP2 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KCP-2; KCP2; Keratinocyte-associated protein 2; KRTCAP2; KTAP2_HUMAN; KTAP2; KRTCAP2 KCP2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Endoplasmic reticulum . Endoplasmic reticulum membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed in skin, heart, placental, liver, skeletal muscle, kidney, pancreas, keratinocytes and dermal fibroblasts.
Function	caution: It is uncertain whether Met-1 or Met-27 is the initiator. . similarity: Belongs to the KRTCAP2 family. . subunit: Component of the oligosaccharyltransferase (OST) complex. OST seems to exist in different forms which contain at least RPN1, RPN2, OST48, DAD1, OSTC, KRTCAP2 and either STT3A or STT3B. OST can form stable complexes with the Sec61 complex or with both the Sec61 and TRAP complexes. . tissue specificity: Expressed in skin, heart, placental, liver, skeletal muscle, kidney, pancreas, keratinocytes and dermal fibroblasts. .
Background	

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products ImagesWestern blot analysis of lysates from PC-12 cells,
primary antibody was diluted at 1:1000, 4°over night