




TM217 rabbit pAb

Catalog No	YP-Ab-08410
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	TMEM217 C6orf128
Immunogen	Synthesized peptide derived from human TM217 AA range: 121-171
Specificity	This antibody detects endogenous levels of TM217 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TM217;TMEM217 C6orf128
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MCF-7 cells,
primary antibody was diluted at 1:1000, 4° over night