




CBY3 rabbit pAb

Catalog No	YP-Ab-08411
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	CBY3
Immunogen	Synthesized peptide derived from human CBY3 AA range: 48-98
Specificity	This antibody detects endogenous levels of CBY3 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CBY3
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	miscellaneous:Chibby is Japanese for 'small'; the gene was so named for the RNAi phenotype seen in flies.,similarity:Belongs to the chibby family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HEK293 cells,
primary antibody was diluted at 1:1000, 4° over night