




RTF1 rabbit pAb

Catalog No	YP-Ab-08418
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	RTF1 KIAA0252
Immunogen	Synthesized peptide derived from human RTF1 AA range: 228-278
Specificity	This antibody detects endogenous levels of RTF1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CDG1N;GTL7;KIAA0252;RNA polymerase-associated protein RTF1 homolog; RTF1;RTF1 KIAA0252
Observed Band	78kD
Calculated Molecular Weight	78kD
Cell Pathway	Nucleus, nucleoplasm .
Tissue Specificity	
Function	function:Plays a role in transcription-coupled histone modification. Plays a role in regulation of transcription. Required for maximal induction of heat-shock genes., similarity:Contains 1 Plus3 domain.,subunit:Interacts with the RNA polymerase II complex. May interact with the PAF1 complex.,
Background	This locus may represent a gene involved in regulation of transcription elongation and chromatin remodeling, based on studies of similar proteins in other organisms. The encoded protein may bind single-stranded DNA. [provided by RefSeq, Sep 2010],
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from HpeG2 cells, primary antibody was diluted at 1:1000, 4° over night