




YN3586

LYPD8 rabbit pAb

Catalog No	YP-Ab-08428
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	LYPD8 UNQ511/PRO1026
Immunogen	Synthesized peptide derived from human LYPD8 AA range: 42-92
Specificity	This antibody detects endogenous levels of LYPD8 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LYPD8;LYPD8 UNQ511/PRO1026
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Lipid-anchor, GPI-anchor . Secreted . Secreted into the lumen of the colon following cleavage of the GPI-anchor. .
Tissue Specificity	Expressed in the large intestine. Preferentially expressed on the epithelial layer exposed to the lumen (at protein level).
Function	similarity:Belongs to the CNF-like-inhibitor family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HpeG2 cells,
primary antibody was diluted at 1:1000, 4°over night