



GABP2 rabbit pAb


Catalog No	YP-Ab-08434
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	GABPB2
Immunogen	Synthesized peptide derived from human GABP2 AA range: 314-364
Specificity	This antibody detects endogenous levels of GABP2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BABPB2;E4TF1;E4TF1-47;E4TF1-53;E4TF1B;GA binding protein transcription factor beta subunit 1 transcript variant gamma 2;GA binding protein transcription factor;beta subunit 1;beta subunit 2;GA-binding protein subunit beta-1;GABP subunit beta-1;GABP subunit beta-2;GABP1_HUMAN;GABPB;GABPB-1;GABPB-2;GABPB1;GABPB2;NRF2B1;NRF2B2;Nuclear respiratory factor 2; Transcription factor E4TF1-47;Transcription factor E4TF1-53;GABP2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	
Function	function:May function as transcription factor capable of interacting with purine rich repeats (GA repeats).,similarity:Contains 5 ANK repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from 293T cells,
primary antibody was diluted at 1:1000, 4° over night