




IQCH rabbit pAb

Catalog No	YP-Ab-08441
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	IQCH
Immunogen	Synthesized peptide derived from human IQCH AA range: 829-879
Specificity	This antibody detects endogenous levels of IQCH at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	IQCH
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Expressed in fetal and adult testis, in germ cells but not somatic cells.
Function	function:May play a regulatory role in spermatogenesis.,similarity:Contains 1 IQ domain.,tissue specificity:Expressed in fetal and adult testis, in germ cells but not somatic cells.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from Hela cells,
primary antibody was diluted at 1:1000, 4° over night