




DEP1B rabbit pAb

Catalog No	YP-Ab-08448
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	DEPDC1B XTP8
Immunogen	Synthesized peptide derived from human DEP1B AA range: 478-528
Specificity	This antibody detects endogenous levels of DEP1B at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DEP1B;DEPDC1B XTP8
Observed Band	
Calculated Molecular Weight	
Cell Pathway	cytosol,
Tissue Specificity	
Function	similarity:Contains 1 DEP domain.,similarity:Contains 1 Rho-GAP domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from A431 cells,
primary antibody was diluted at 1:1000, 4° over night