




COX16 rabbit pAb

Catalog No	YP-Ab-08452
Isotype	IgG
Reactivity	Human; Mouse
Applications	IHC;IF
Gene Name	COX16 C14orf112 HSPC203 PTD019
Immunogen	Synthesized peptide derived from human COX16 AA range: 9-59
Specificity	This antibody detects endogenous levels of COX16 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:50-200. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C14orf112;COX 16;hCOX16;HSPC203;MC4DN22;COX16;COX16 C14orf112 HSPC203 PTD019
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Mitochondrion inner membrane ; Single-pass membrane protein .
Tissue Specificity	Widely expressed. Expressed at higher level in skeletal muscle, heart and liver.
Function	miscellaneous:No COX16 mutations have been detected in patients with cytochrome c oxidase (COX) deficiency.,similarity:Belongs to the COX16 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 30min).