




GRIFN rabbit pAb

Catalog No	YP-Ab-08455
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	GRIFIN
Immunogen	Synthesized peptide derived from human GRIFN AA range: 75-125
Specificity	This antibody detects endogenous levels of GRIFN at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GRIFN; GRIFIN
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Not detected in lens.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MCF-7 cells,
primary antibody was diluted at 1:1000, 4° over night