




5NT3L rabbit pAb

Catalog No	YP-Ab-08456
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	NT5C3L
Immunogen	Synthesized peptide derived from human 5NT3L AA range: 23-73
Specificity	This antibody detects endogenous levels of 5NT3L at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	5NT3L;NT5C3L
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm .
Tissue Specificity	
Function	catalytic activity:A 5'-ribonucleotide + H(2)O = a ribonucleoside + phosphate., function:Can act both as nucleotidase and as phosphotransferase., similarity:Belongs to the pyrimidine 5'-nucleotidase family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from Hela cells,
primary antibody was diluted at 1:1000, 4° over night