



K1614 rabbit pAb


Catalog No	YP-Ab-08461
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	KIAA1614
Immunogen	Synthesized peptide derived from human K1614 AA range: 775-825
Specificity	This antibody detects endogenous levels of K1614 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	63441H18Rik;Coiled coil domain containing protein MTMR15;DKFZp451H236; DKFZp686K16147;FAN1;FAN1_HUMAN;FANCD2 and FANCI associated nuclease 1;FANCD2/FANCI associated nuclease 1;FANCD2/FANCI-associated nuclease 1;Fanconi associated nuclease 1;Fanconi-associated nuclease 1; hFAN1;KIAA118;KMIN;MTMR 15;MTMR15;Myotubularin related protein 15; Myotubularin-related protein 15;K1614;KIAA1614
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from CACO2 cells,
primary antibody was diluted at 1:1000, 4° over night