




LUZP2 rabbit pAb

Catalog No	YP-Ab-08466
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	LUZP2 UNQ2566/PRO6246
Immunogen	Synthesized peptide derived from human LUZP2 AA range: 135-185
Specificity	This antibody detects endogenous levels of LUZP2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LUZP2;LUZP2 UNQ2566/PRO6246
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	
Function	disease:Deletion of the LUZP2 gene may be a cause of Wilms tumor--aniridia--genitourinary anomalies--mental retardation syndrome (WAGR) [MIM:194072].,
Background	This gene encodes a leucine zipper protein. This protein is deleted in some patients with Wilms tumor-Aniridia-Genitourinary anomalies-mental Retardation (WAGR) syndrome. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Oct 2011],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HpeG2 cells,
primary antibody was diluted at 1:1000, 4°over night